

CRYSTAL M100

DENTAL UNIT



MEDTECH
REMEDIES PVT. LTD.

Ultimate Comfort

A stylish and smart dental chair that will infuse comfort and confidence in you as well your patients. This ergonomically designed chair will put you at ease

Key Features

- Smart designed dental chair with 24v noiseless DC motor, with 3 memory programmes
- Wide Instrument working tray
- Under hanging delivery unit with 5 point inclusive of one micromotor point
- Movable hydro group for easy maintenance.
- Large ergonomic soft PU cushion
- Reflective LED light with intensity of 32k Lux controlled with 6 stages
- Ceramic Spittoon with 90 degree movement
- Auto spittoon flush and cup filler control system
- Hot and cold water control system
- LED film Viewer
- Multifunctional foot control with Chip Air/Water Activation
- Saliva ejector with high and low pneumatic suction



Major Components highlights




Rotatable Spittoon Unit



Ergonomic Design PU Upholstery



Assistant Instrument Table



New Design LED Lamp



Dentist Stool

CRYSTAL M100

Standard Unit

- Large Ergonomic Soft PU Cushion
- Motorised dental chair control system - 24V Noiseless DC under voltage motor chair
- Reflective LED operating light
- Movable ceramic spittoon
- Water supply system
- DC under voltage X-ray Film Viewer
- Programmable cup filler and bowl rinsing
- Water heater
- Air turbine hose without light 2 sets
- Air motor hose without light 1 set

Upgrades

- Scaler and Light cure
- Brushless Micromotor
- Intraoral Camera And 17' HD Monitor
- 1 Pc Airmotor
- 2 Pcs Air Turbine
- Leather Cushion With memory foam

Technical Parameters

- Power Supply: 220V, 50/60Hz
- Air Source Pressure: 0.5-0.8MPa
- Water Source Pressure: 0.2-0.8MPa
- Motor Voltage : 24V



Standard Upholstery Colors





MEDTECH
REMEDIES PVT. LTD.



+91 8590024281



india@med-tech.com



med-tech.in



No C-3, 2nd Floor,
KSSIDC Industrial Estate, Yelahanka,
Bangalore - 560064, Karnataka, India.



f @ medtechremedies

Dental care through innovation



PLANMECA®



mectron

ALERIO®
X-RAYS

COXO®