

MED AIR

COMPRESSOR UNIT



MEDTECH
REMEDIES PVT. LTD.

Dental Compressor Unit

Medtech's Oil Free Compressors deliver high performance and consistently deliver clean, dry, oil-free, hygienic air. This compressor consumes low energy and is fully automatic, automatic start and stop control. Our dental air compressors are robust, powerful, durable, and efficient. Clean dry air not only protects your patients and the work you perform for them, but also your valuable equipment as well reduced maintenance and longevity of the unit is increased.

MED AIR 100

Specifications

- Voltage / Frequency : 220/230 V 50/60 Hz
- Free air delivery : 135 L/min
- Tank capacity : 30 L / 40 L
- ON/OFF pressure : 8 bar
- Power : 1 HP
- Noise level : 60dB(A) at 1m distance



MED AIR 150

Specifications

- Voltage / Frequency : 220/230 V 50/60 Hz
- Free air delivery : 200 L/min
- Tank capacity : 50 L
- ON/OFF pressure : 8 bar
- Power : 1.5 HP
- Noise level : 70dB(A) at 1 m distance

For more details:

+91 8590024281

india@med-tech.com

med-tech.in

@medtechremedies