

MED VAC

SUCTION UNIT



MEDTECH
REMEDIES PVT. LTD.

Dental Suction Unit

A dental aspiration system is one of the essential elements of any clinic or dental practice. Its main use is related to the reduction or elimination of cross contamination between the patient and the dentist as well as maintaining a hygienic environment for the clinic. The MedVac suction unit is a powerful and effective aspiration with air/fluid separation system that is made for continuous use; thereby giving a comfortable and less time consuming dental procedures for both the Dentist and the patient.

MED VAC

Specifications

- Model : MED VAC (With Collection Jar)
- Type: Mobiles
- Voltage / Frequency : 220/230 V 50/60 Hz
- Air flow: 85 (m³/h)
- Pressure : 120 (mbar)
- Vacuum : -110 (mbar)
- Power : 0.25 KW
- Noise level : <55 dB(A) at 1 m distance



For more details:

+91 8590024281

india@med-tech.com

med-tech.in

@medtechremedies