

CLASSIC M50

DENTAL UNIT



Best Treatment Experience

A dental unit that proclaims safety, Comfort and durability. With good features and easy to use controls, you will never stress yourself and your patients

Key Features

- Electrically operated dental chair with 24v noiseless DC motor, with 2 memory programmes
- Under hanging delivery unit with 5 point inclusive of one micromotor point
- Halogen Lamp or Sensor based LED light with 2 stages
- Ceramic Spittoon with 90 degree movement, Cup filter control system
- Cup filler with both hot & cold water control system
- LED film Viewer
- Multifunctional foot control
- Saliva ejector with high and low pneumatic suction



Major Components highlights



PU leather Cushion



Dentist Stool



Ceramic Cuspidor



Halogen / LED Operating Lamp



Concealed PC Box



Assistant Tray

Standard Unit

- Underhanging type of Instrument version
- Polyether air and water tubing
- Integrated Multifunctional foot control
- 90° Rotatable ceramic spitton
- Multifunctional assistant control box
- Built out bottle for water supply or disinfectant
- Water heater
- 24V Noiseless DC under Voltage motor chair
- Halogen / LED Operating lamp
- 3 way Syringe

Upgrades

- Leather cushion
- High speed air turbine hand piece
- Low speed air turbine hand piece
- Intra Oral Camera
- Ultrasonic scaler
- Curing light
- Inbuilt micromotor

Technical Parameters

- Power Supply: 220V, 50/60Hz
- Air Source Pressure: 0.55-0.6MPa
- Water Source Pressure: 0.2-0.4MPa
- Motor Voltage : 24V



Standard Upholstery Colors





MEDTECH
REMEDIES PVT. LTD.



+91 8590024281



india@med-tech.com



med-tech.in



No C-3, 2nd Floor,
KSSIDC Industrial Estate, Yelahanka,
Bangalore - 560064, Karnataka, India.



f @ in @medtechremedies

Dental care through innovation



PLANMECA[®]



mectron

ALERIO[®]
X-RAYS

COXO[®]