



Image Plate Scanner **EQ600**

Professional & High Quality Image

VRN Image Plate Scanner Specification			
Data	Size	Standard	High
Support IP sizes	0 (22×31mm)	343×484 pixels	628×885 pixels
	1 (24×40mm)	375×625 pixels	685×1143 pixels
	2 (31×41mm)	484×640 pixels	886×1171 pixels
	3 (27×54mm)	421×843 pixels	771×1542 pixels
Product Brand			
Product Model	EQ-600		
Selectable Pixel Size	35µm		
Detector Technology	APS CMOS		
Software*	UT-Dental		
Active Area	32*54mm		
Mechanical Size	34*55mm		
Readout Time	15s		
Data Capture	16-Bits 65,536 grey level		
Eraser	Embedded(auto)		
Dimensions (H*W*D)	305*211*182mm		
Weight	6Kg		
System Configuration	Windows XP,Windows 7,Windows 8,Windows 10		
Power	100-240V(AC) 50/60Hz 15A		
Interface	USB 2.0 high speed(480 Mbps) / Ethernet(100Mbps)		
Work Temperature	10 to 40 °C		
Storage Temperature	-25 to +60°C		
Humidity (non-condensing)	10 to 93% R.H		

