

Scaler & Lightcure

In-built Scaler | Table top | Lightcure



Technical Specifications

- Description: In-built scaler
- Model : I 02 / I 02E
- Detachable handpiece
- I 02 : Scalling & Perio
- I 02E : Scalling & Perio & Endo
- specifications same as K08A
- Dim: 55mm x 160mm x 160mm



Technical Specifications

- Description: Table top scaler
- Model : K08A
- Detachable Handpiece
- Scaling & Perio
- Tips : G1+2, G2, G4, P1
- Water pressure: 0.01Mpa- 0.05Mp
- Vibration frequency of the tip: 25-31 KHz
- Ultrasonic output power: 3w-20w
- Input power : 30-48 VA
- Main unit input: DC 30V 1.5A



Technical Specifications

- Description: Lightcure
- Model: V200
- Light intensity: 900 - 1650W/cm
- Device type: cordless exchangeable Li-ion battery
- Working time : 3s,5s,10s,15s,20s,25s,30s,35s,40s
- Working mode : full power, progressive, pulse



For more details:

+91 8590024281

india@med-tech.com

med-tech.in

@medtechremedies